

PHOTON SOUND BEAM

Rife-Lakhovksy-Lymphatic

PlasmaElectrical-RF Product Line



Photon Sound Beam Infinity RF

Clinical Model

Electrical & RF Plasma Tubes, Dual-Oscillating RF Sound Probe

Why the Photon Sound Beam is unique

The *Photon Sound Beam* is unique in the rapidly emerging arena of Electro-medicine. The *Photon Sound Beam* is actually several technologies in one unit – utilizing light, sound, subtle energy, plasma-electricity and radio-wave delivery. This “Vibrational Energy Medicine” understands the simple fundamental nature of health and wellness as a vibrational relationship with the intelligence of the body. Frequency is the language of energy, all matter vibrates as energy; indeed, matter is nothing but energy. More specifically, matter is condensed light carried into physical manifestation by photons that are directed by creative intelligence.

The focus of Plasma Electrical Technology like the *Photon Sound Beam* is to provide naturalness, wholeness and cleansing with the broadcast of pure-waveform frequencies and frequency harmonics. In this way, the body is better prepared to receive frequency information required for the ideal resonance and energy state of each cell.

Fundamentally “health” is simply a matter of what gets into the cell and what is eliminated from the cell. Electro-chemistry plays a key role here. It orchestrates the transport of nutrients and the elimination of toxins – and this is directly related to the electrical properties of the cell.

Plasma-electricity is natural to the body. Plasma-electrical Instruments deliver natural biological current deep into the body at the level of the cell. It quickly re-establishes the ideal electrical energy state and frequency of every cell in the body. It removes blockage and barrier issues by “electrical disassociation”, opening-up the body’s energy/information and elimination channels. Plasma-electricity is a perfect way to deliver Rife and other frequencies in a way which is bio-available.

Plasma-Electricity from the *Photon Sound Beam* results in:

- **Deep cleansing and oxygenation of tissues/fluids** enhancing the flow of the Lymphatic and Cardiovascular Systems
- **Removal of "Blockage & Barrier Issues"** created from acid-forming diets while enhancing the flow of all the natural elimination channels
- **Improved delivery & integration (bio-availability) of frequencies** to devitalize pathogens & parasites while improving the health of normal cells in the body
- **Photobiotic nourishment** from light and subtle energies to enhance the Endocrine System

“Our physiology is primarily electrical in nature, and secondarily physical.”

This is why we specialize in Plasma-Electrical Instruments like the ***Photon Sound Beam***. The ***Photon Sound Beam*** uses hand-held **Electric Plasma Tubes** -- glass tubes that are ionized by Tesla Coils and filled with Noble Gases (Xenon/Argon/Krypton) in a near vacuum state. Plasma-electricity deeply penetrates the body while cleansing and oxygenating all the tissues and fluids while making the delivery of frequencies more bio-available for the body. Note: See Electric Plasma Tubes below.

At the same time the ***Photon Sound Beam***’s ***RF Sound Probe*** component penetrates the frequencies with high-harmonic content to re-establish the ideal resonance of the cell. These RF Carriers contain the same frequencies emitted by the RF-enhanced Electric Plasma Tubes is typically used synergistically at the same time, however one may use ***RF Sound Probe*** by itself especially for night time with close distance (1-3 ft.) broadcasting. This component has been vastly upgraded and uses a highly scalar braided, spring-shaped broadcast antenna to maximize its subtle energy and its overall holistic value.

Finally, the ***Photon Sound Beam*** has advanced control features such as Frequency Selection and Individual Pause & Pulse -- to overcome body accommodation and to insure long-term benefits.

From our longtime personal experience and hundreds of ***Photon Sound Beam*** users since 1994 - this evolved combination of many technologies *into one easy-to-use instrument* comprehends and out-performs any technology available today, including: Rife RF-Broadcast or Electric Metal Rod/Pad Instruments, Lakhovsky Multiwave Oscillators (MWO’s), Pulsed EMF, Light Therapy, any hand-held instrument or any other Plasma Tube Instrument. Unlike several other manufacturers -- We offer lifetime technical support and a full 2-year warranty.



Photon Sound Beam XII
Most Popular Personal-Use & Clinical Model

Plasma-Electrical Lymphatic Rife/Lakhovsky Technology

Since 1988 we have been educating, creating and marketing a wide variety of the world's best Vibrational Energy Instruments. We have been evolving *Plasma-Electrical/RF Technology* since Ed Skilling's original 1994 single frequency *Photon Sound Beam*.

The ***Photon Sound Beam*** is a *unique wrap-around technology* combining the benefits of many instruments all into one unit. It is called a "**Plasma-Electrical Lymphatic Rife/Lakhovsky Technology**" -- it combines the remarkable *deep tissue-cleansing benefits* of Lymphatic Drainage Technology with more aggressive and more allopathic Rife & Lakhovsky Technologies which use frequencies and frequency harmonics to devitalize harmful pathogens and parasites.

The NEW ***Photon Sound Beam*** has a vastly upgraded circuitry and construction with EXPANDED CONTROL FEATURES such as Frequency Selection and Individual Pause & Pulse Control. Pulsation Control is essential for enhancing regeneration in the body. In general, these expanded control features help in overcoming "body accommodation" which happens when the body's intelligence begins to ignore repeated frequency information and energy.

Expanded control features can meet the needs of a variety of individuals and conditions. It is a key factor in the fast-acting universal response by *Photon Sound Beam* users.

The *Photon Sound Beam* now includes over 20 upgrades and comes in 4 different models.

All *Photon Sound Beam* models include:

- **Plasma-Tube Glass Applicators** with Xenon, Argon, and Krypton Noble Gases for removal of blockages and barrier issues
- **Unique RF Sound Probe** component with its high-harmonic RF Carriers & Rife Frequencies generated from a dual-oscillating quartz crystal oscillators and based on the work of Georges Lakhovsky

[PSB Model Comparison](#)

Electric Plasma Tubes & Lymphatic/Tissue Cleansing

It is well known that ONLY *Electrical Plasma Tubes* – not *Radio wave Plasma Tubes* – can provide *deep cellular disassociation of all blockages and barrier issues*. This "Plasma-Electricity" deeply penetrates the body and follows its natural electrical circuitry. However this extremely low current *cannot* travel along electrical wires. It is airborne, electrostatic electricity and it must be contained inside the ionized Noble Gas atmosphere of glass tubes applied directly to the body. Electric Plasma Tubes provide an excellent media for transmitting frequencies and they deeply cleanse the lymph, fluids and tissues. They also provide: (1) the elemental colors frequencies of the ionized Noble Gases (Xenon/Argon/Krypton) and (2) highly scalar Zero-point Subtle Energies. The plasma-electricity emitted from Electrical Plasma Tubes functions at biological levels of electrical current, which is gentle and diffuse rather than laser-like as with metal mod delivery. It automatically dis-associates bound cells, bound protein and toxins while encouraging the flow of all the natural elimination channels in the body and deeply cleansing all the tissues and fluids.

This cannot occur with either type of popular Rife Technologies used today: Rife Metal Rod Technology or Rife Broadcast Technology. Metal Rod output cannot effectively delivery frequencies because of the high skin impedance. Rife Broadcast instruments emit frequencies vulnerable to body accommodation and environmental EMF pollution. Neither instrument effectively provides any lymphatic drainage or enhanced cellular integration of frequencies found with the *Photon Sound Beam*. NOTE: Both types of Plasma Tube Output are included in our clinical *PSB Infinity RF*, which also has 60 frequency programs and an additional OxyGold Triple Gas Plasma Tube applicator.

Some instruments such as the *Photon Genie* -- a single frequency knock-off of the 1994 Prototype *Photon Sound Beam* -- use off-the-shelf painted glass applicators used in Neon signs, which block the transmission. Our *Photon Sound Beam* incorporates 3 Noble Gases in a specific proportion, for a pure transmission of the elemental colors of Xenon, Argon and Krypton and their exotic healing value. Noble Gases, ionized in vacuum with biological levels of electricity, is an ideal media to emit frequencies, subtle energy and *photobiotic nourishment*. It is the same media found in intergalactic space and it is the basis of the quantum value of our *Photon Sound Beam* PlasmaElectrical-RF Product Line.

[More PSB Information](#)

Plasma-Electrical Instruments: A Holistic Approach to Energy Medicine

The *Photon Sound Beam* is a holistic Vibrational Energy Technology – it offers the best of both worlds. We have found that in a toxic world -- effective long-term benefits with any instrument must involve a balance between allopathic and holistic approaches to health. When delivering energy and information to the body -- NOURISHMENT MUST BE A PART OF ANY APPROACH TO HEALTH for faster acting and long-term results.

Electro-medical Instruments which deliver energy and frequencies to the body should have the same considerations as food – wholeness, naturalness and freshness:

Naturalness -- they should deliver natural biological levels of electrical current. One of the most important factors in creating an effective holistic Vibrational Energy Technology is to use an *energy carrier*, which can effectively travel along the body's natural circuitry to reach every cell.

Wholeness -- the frequencies should have high-harmonic content. In the 1920's Georges Lakhovsky's overall conclusion using his *MultiWave Oscillator* experimenting on cancer cells in plants and animals was -- that it did *not* matter which Radio wave (Rife) Frequency he used -- but that the *frequency with the highest amount of harmonics* always worked the best. In other words, using the exact "frequency per dis-ease" is not as important as using a bio-interactive frequency *rich in harmonic content*. In addition the frequencies must be effectively delivered and received by the body.

Freshness -- they should have control features such as variable pulsation and frequency selection to avoid accommodation by the body. A truly effective instrument must *factor* the ability of the body's intelligence to accommodate or ignore any broadcast of frequencies and energy -- an expression of protective homeostasis.

Fundamentally health is simply a matter of what gets into the cell and what is eliminated from the cell. Electro-chemistry plays a role here. The "software of the cell" is its outer membrane, which orchestrates the transport of nutrients and the elimination of toxins -- and this is directly related to the electrical properties of the cell.

It is well established that acid-forming diets (low pH levels) contribute to the clustering of red blood cells and other cells as well as congealed lymph. Every function of our bodies relies on these basic electrical properties. And every issue of health is directly tied to cellular electrical energy states.

The *Photon Sound Beam's* Plasma Tube applicators function on the level of biological levels of electricity (extremely low micro-current) which quickly re-establishes an "ideal cellular resonance" and an "ideal electrical polarity". This is the result of a process called ELECTRICAL DIS-ASSOCIATION. Dark Field Microscopy (Live Cell Analysis) has repeatedly shown an 80% improvement in the health of red blood cells in one 20 minute session with the *Photon Sound Beam*.

ELECTRICAL DIS-ASSOCIATION: the separation of bound protein matter, bound cells and toxins by the action of biological micro-current on the cell. The cell membrane has an electrically negative (-) charge and the cell nucleus is positively (+) charged. Electrical energy is a flow of negatively (-) charged electrons which can quickly re-establishes a (-) polarity to this vital membrane software of the cell.

Evolution of the *Photon Sound Beam*

Our *Photon Sound Beam* Product Line combines the best of early pioneers and researchers who made major contributions to 20th Century Electromedicine.

In the early part of the 20th Century researchers experimented with Radio wave and Electrical Fields and discovered dramatic effects in treating serious issues of health.

In the 1920's, after 20 years of research experimenting with Radio waves to treat cancer in plants and humans, Georges Lakhovsky discovered that Radio wave Frequencies with high harmonic content exhibited remarkable success in short period of time. He is famous for his Multiwave Oscillator, which broadcasted powerful electric fields at distance.

Later in the 1930's -- Royal Raymond Rife identified and demonstrated specific "mortal oscillatory rates" or frequencies associated with pathogens responsible for cancer, tuberculosis and other life-threatening diseases. With the help of Lee DeForest (the "father of modern electronics") Rife then developed a Radio wave Plasma "*Beam Tube*" which was found to quickly devitalize these pathogens at a distance with a powerful shortwave radio wave broadcast carrying mortal oscillatory rates.

In the late 1970's our current business partner introduced an electrical "Metal Rod Rife Frequency Generator" (created by John Crane, Rife's assistant in the 1950's) to Ed Skilling and Robert Beck. This Electrical Frequency Generator was advanced by in the 1980's and experimented in clinics in Mexico. Since the specifications were made public --today there are many vendors of Metal Rod Rife Frequency Generator with many frequency applications.

However -- when we began working with Ed Skilling in the late 1980's -- he mentioned that instruments, which use Metal Rods, deliver a *higher current* than what is natural to the body and can only follow the main conductive pathways of the cardio-vascular system. Skilling asserted that this type of Rife Technology has limited penetration and *distorts the waveforms of specific Rife Frequencies* due to metal-to-skin "impedance". This prompted Ed to create his *Sound Probe* in 1988, which broadcasts a radio wave signal with RF Quartz Crystal Oscillators. The *Sound Probe* is applied close to the body to avoid chaotic EMF and cell phone pollution found with Rife broadcasting at distance.

So In 1994, Ed Skilling developed his "*Inert Gas Sound Probe*" which added Glass Applicators filled with Inert Gases. These Electrically Ionized Plasma Tubes have a significant advantage over Rife's *radio wave* Plasma Tubes because they can deeply penetrate *biological levels of Plasma-electricity* that deeply "cleanse" the Lymphatic System, Tissues and other fluids such as cerebral spinal fluid and cellular fluids. Working with Ed Skilling we re-named this instrument the *Photon Sound Beam*.

In 1997 our very special and knowledgeable engineers began to evolve this 94 prototype *Photon Sound Beam* by creating the *PSB VII, VII and Essential PSB*. *These advancements included:*

(1) dramatically cleaning-up the circuitry, (2) improving the durability of the hardware, and (3) adding important basic control features such as Universal Rife Frequency Selection *and* Pulsation Control. These enhanced control features are important for enhancing regeneration in body, and for overcoming a common problem with any technology: *body accommodation* --

when after time the body ignores any energy or frequencies. This typically happens after 30 to 40 minutes exposure to any Vibrational Energy Technology.

In 2004 we made a quantum advancement of the *Photon Sound Beam* technology by adding a second Tesla Coil so that both Electric Plasma Tubes output electrical energy rather the previous Mono-polar *Photon Sound Beams*. This allows an omni-directional flow of biological micro-current from the Noble Gas glass applicators as well as enhancing cleansing action and the frequencies they carry. In addition, we added sophisticated computer chip circuitry. This enabled frequency ranges from 1 to 10,000 hertz, which covers most of the Rife Frequency Range.

TODAY the [*Photon Sound Beam XII*](#) is our most popular model. It is both a personal-use and clinical instrument; it has 12 Brainwave and Rife Frequency Selections and a 36 Priority Rife Frequency Sweep, Individual Pause and Pulse Controls, and vastly improved reliability and appearance.

Our clinical model -- [*Photon Sound Beam Infinity RF*](#) -- adds full-frequency programming with up to 60 Frequency Programs including 40 Pre-Set Programs. The *PSB Infinity* includes a (1) Rife Frequency book ([The Electroherbalism Frequency Lists](#)), an additional glass applicator and BOTH Electrical and Radiowave Plasma Tube output.

PHOTON SOUND BEAM

Product Line

You can trade up to any model!

All *Photon Sound Beam* models come with 2-3 glass Plasma Tubes & the RF Sound Probe. The RF Sound Probe is the 2nd component of the *Photon Sound Beam*. It emits 2 harmonically-rich RF Signals in Phi Ratio (ubiquitous in Nature) to carry the same frequencies coming out the Glass Tubes. The RF Sound Probe attachment can be used either separately or at the same time as the electrical glass Plasma Tubes.

Essential Photon Sound Beam – This is a much upgraded version of the original 1994 *PhotonSoundBeam* prototype still being manufactured by other sources (such as the *Photon Genie* at twice the price). It has a *Monopolar* Plasma Tube Output, Pulse Rate Control and one pre-set Primary Universal Rife Frequency (728 Hz) rich in bio-available harmonics. A Monopolar System means that one glass Plasma Tube emits energy and the other one acts as a ground to complete the circuit. \$1495

Photon Sound Beam XII -- This is our most popular personal use/clinical model and utilizes 2 Tesla Coils resulting in a Bi-Phasic Output. This means that *both* glass applicators emit energy and frequencies – resulting in an omni-directional output. Each applicator RF Enhanced and is independent of the other, and the biphasic output *eliminates any annoying electrostatic shocking*. The *PhotonSoundBeam XII* has a 36 Sweep Setting that automatically cycles all of the major Rife and Brainwave frequencies *each pulsed 5 times*. It also has 4 Universal Rife Frequencies, 5 classical cancer frequencies and 2 brainwave frequency settings as well as individual Pause and Pulse Controls and a vastly upgraded durability and appearance. **\$2395**

PSB Infinity RF –This clinical model has 3 components – three instruments in one. It combines the *Electrical PlasmaTube Ionizer* and the *RF Sound Probe* of the *PhotonSoundBeam XII* with a *Radiowave PlasmaTube Ionizer*. The *PSB Infinity RF* has full selection of frequencies -- 60 Frequency Programs containing up to 20 frequencies each and 40 Preset Programs. It comes with an additional 3rd (OxyGold) Glass Applicator with Oxygen and Gold resonance. The application of BOTH the ELECTRICAL and RADIOWAVE IONIZATION provides deep tissue/fluid cleansing and oxygenation – as well as a powerful synergy with intimate, nourishing Rife Broadcast Technology. The Radio wave Plasma Tubes of the *PSB Infinity RF* deeply penetrate frequencies while the Electrical Plasma Tubes & RF Sound Probe deeply cleanse and nourish the cells and tissues – while continuing to address pathogens and at a deeper level. **\$3595**

The ***LifePulse Pro*** is the same as the *Infinity RF* except that it ADDS Electrical Metal Rod/ Plate Rife Technology. *This model is a special order item and has limited availability.* **\$4795**

Please call for a free introductory consultation, or to answer any questions.

Linda and Keith -- Integrated Healing Systems -- Unified Field Technologies
www.PhotonSoundBeam.net – healthyvibes@cableone.net 928 634 0833

Important Medical Disclaimer: The Photon Sound Beam is for experimental and research purposes only. No medical claims are expressed or implied.